

Meas. No. 4153 SOP 0 Date 04.08.2008 15:12:00 Operator FRITSCHLAN\Benes

Material: **PVC Powder**

L.V.M.

Description: Settings for dry measurement: pressure 3bar, beam obscuration : 2-4%

Mode	Dry	Serial No.	22.2000.00/90771
Beam absorption	2.8 %	Cellposition	1
Meas. range	4.93 µm - 677.07 µm	Calculation	Automatische Modellerkennung
Channels	51	Prot.-No.:	080262
Scans	100		

5-99%

Obere Kornklasse [µm]	Q3(x) [%]
100.264	5.0
110.477	10.0
117.176	15.0
123.212	20.0
127.981	25.0
132.647	30.0
137.224	35.0
141.535	40.0
145.846	45.0
150.183	50.0
154.909	55.0
159.636	60.0
164.362	65.0
170.278	70.0
176.382	75.0
182.907	80.0
192.134	85.0
202.404	90.0
218.694	95.0
248.006	99.0

01-600µm

Obere Kornklasse [µm]	Q3(x) [%]
0.150	0.0
0.500	0.0
1.000	0.0
3.000	0.0
5.000	0.0
8.000	0.1
10.000	0.1
15.000	0.1
20.000	0.2
30.000	0.2
40.000	0.3
50.000	0.6
60.000	0.7
80.000	1.4
100.000	4.9
150.000	49.8
200.000	89.3
250.000	99.1
300.000	100.0
400.000	100.0
500.000	100.0
600.000	100.0

