

Meas. No. 6790 SOP 45 Date 05.05.2009 14:19:13 Operator: FRITSCHLAN\Benes

Material: limestone

P7pl 80ml 1mm 30min

Description: water + 30s ultrasonic + 0.1% Na4P2O7

Calculation Automatische Modellerkennung

Mode Wet

Serial No. 22.2000.00/90771

Beam absorption 13,0 %

Pump 70,00 %

Meas. range 0,10 µm - 53,10 µm

Cellposition 1

Channels 51

Ultrasonic On

Scans 100

Prot.-No.: **090082**

5-99%

Obere Kornklasse [µm]	Q3(x) [%]
0,172	5,0
0,209	10,0
0,240	15,0
0,272	20,0
0,302	25,0
0,334	30,0
0,367	35,0
0,402	40,0
0,438	45,0
0,478	50,0
0,521	55,0
0,567	60,0
0,621	65,0
0,682	70,0
0,752	75,0
0,838	80,0
0,947	85,0
1,099	90,0
1,342	95,0
1,874	99,0

001-5 µm

Obere Kornklasse [µm]	Q3(x) [%]
0,010	0,8
0,050	0,8
0,100	0,8
0,200	8,7
0,300	24,7
0,400	39,7
0,500	52,7
0,600	63,1
0,700	71,4
0,800	78,0
0,900	83,1
1,000	87,0
1,500	96,9
2,000	99,3
2,500	99,9
3,000	100,0
3,500	100,0
4,000	100,0
4,500	100,0
5,000	100,0

