

Meas. No. 6841 SOP 45 Date 15.05.2009 09:30:46 Operator: FRITSCHLAN\Benes

Material: limestone

P7pl 80ml 1mm 60min

Description: water + 30s ultrasonic + 0.1% Na4P2O7

Calculation Automatische Modellerkennung

Mode Wet

Serial No. 22.2000.00/90771

Beam absorption 11,8 %

Pump 70,00 %

Meas. range 0,10 µm - 53,10 µm

Cellposition 1

Channels 51

Ultrasonic On

Scans 100

Prot.-No.: **090082**

5-99%

Obere Kornklasse [µm]	Q3(x) [%]
0,187	5,0
0,223	10,0
0,254	15,0
0,283	20,0
0,311	25,0
0,340	30,0
0,370	35,0
0,400	40,0
0,432	45,0
0,466	50,0
0,503	55,0
0,543	60,0
0,587	65,0
0,636	70,0
0,695	75,0
0,764	80,0
0,851	85,0
0,969	90,0
1,156	95,0
1,549	99,0

01-10µm

Obere Kornklasse [µm]	Q3(x) [%]
0,100	0,4
0,200	6,7
0,300	23,1
0,400	39,9
0,500	54,7
0,700	75,5
1,000	91,1
1,500	98,8
2,000	99,9
2,500	100,0
3,000	100,0
4,000	100,0
5,000	100,0
6,000	100,0
7,000	100,0
8,000	100,0
9,000	100,0
10,000	100,0

