

Meas. No. 6782 SOP 45 Date 05.05.2009 08:03:32 Operator: FRITSCHLAN\Benes

Material: limestone

P7pl 80ml 2mm 120min

Description: water + 30s ultrasonic + 0.1% Na4P2O7

Calculation Automatische Modellerkennung

Mode Wet

Serial No. 22.2000.00/90771

Beam absorption 7,9 %

Pump 70,00 %

Meas. range 0,10 µm - 53,10 µm

Cellposition 1

Channels 51

Ultrasonic On

Scans 100

Prot.-No.: **090082**

5-99%

Obere Kornklasse [µm]	Q3(x) [%]
0,209	5,0
0,248	10,0
0,284	15,0
0,316	20,0
0,350	25,0
0,383	30,0
0,417	35,0
0,452	40,0
0,490	45,0
0,530	50,0
0,573	55,0
0,620	60,0
0,674	65,0
0,733	70,0
0,802	75,0
0,887	80,0
0,994	85,0
1,140	90,0
1,381	95,0
1,886	99,0

001-5 µm

Obere Kornklasse [µm]	Q3(x) [%]
0,010	0,1
0,050	0,1
0,100	0,1
0,200	4,2
0,300	17,5
0,400	32,5
0,500	46,3
0,600	57,9
0,700	67,4
0,800	74,9
0,900	80,7
1,000	85,2
1,500	96,6
2,000	99,3
2,500	99,9
3,000	100,0
3,500	100,0
4,000	100,0
4,500	100,0
5,000	100,0

