

Meas. No. 6767 SOP 45 Date 04.05.2009 10:21:23 Operator: FRITSCHLAN\Benes

Material: limestone

P7pl 80ml 2mm 30min

Description: water + 30s ultrasonic + 0.1% Na4P2O7

Calculation Automatische Modellerkennung

Mode Wet

Serial No. 22.2000.00/90771

Beam absorption 12,5 %

Pump 70,00 %

Meas. range 0,10 μm - 53,10 μm

Cellposition 1

Channels 51

Ultrasonic On

Scans 100

Prot.-No.: **090082**

5-99%

Obere Kornklasse [μm]	Q3(x) [%]
0,254	5,0
0,307	10,0
0,353	15,0
0,396	20,0
0,437	25,0
0,479	30,0
0,522	35,0
0,565	40,0
0,612	45,0
0,662	50,0
0,714	55,0
0,774	60,0
0,840	65,0
0,911	70,0
0,999	75,0
1,103	80,0
1,233	85,0
1,410	90,0
1,689	95,0
2,307	99,0

01-10 μm

Obere Kornklasse [μm]	Q3(x) [%]
0,100	0,1
0,200	1,5
0,300	9,2
0,400	20,5
0,500	32,5
0,700	53,6
1,000	75,0
1,500	92,1
2,000	97,7
2,500	99,5
3,000	99,9
4,000	100,0
5,000	100,0
6,000	100,0
7,000	100,0
8,000	100,0
9,000	100,0
10,000	100,0

