

Meas. No. 6778 SOP 45 Date 04.05.2009 15:43:15 Operator: FRITSCHLAN\Benes

Material: limestone

P7pl 80ml 2mm 60min

Description: water + 30s ultrasonic + 0.1% Na4P2O7

Calculation Automatische Modellerkennung

Mode Wet

Serial No. 22.2000.00/90771

Beam absorption 14,0 %

Pump 70,00 %

Meas. range 0,10 µm - 53,10 µm

Cellposition 1

Channels 51

Ultrasonic On

Scans 100

Prot.-No.: **090082**

5-99%

Obere Kornklasse [µm]	Q3(x) [%]
0,236	5,0
0,283	10,0
0,323	15,0
0,361	20,0
0,398	25,0
0,435	30,0
0,473	35,0
0,513	40,0
0,553	45,0
0,598	50,0
0,645	55,0
0,697	60,0
0,755	65,0
0,819	70,0
0,896	75,0
0,989	80,0
1,104	85,0
1,261	90,0
1,509	95,0
2,052	99,0

01-10µm

Obere Kornklasse [µm]	Q3(x) [%]
0,100	0,2
0,200	2,2
0,300	12,1
0,400	25,2
0,500	38,5
0,700	60,3
1,000	80,6
1,500	94,9
2,000	98,8
2,500	99,8
3,000	100,0
4,000	100,0
5,000	100,0
6,000	100,0
7,000	100,0
8,000	100,0
9,000	100,0
10,000	100,0

