

Meas. No. 7092 SOP 18 Date 18.06.2009 10:01:05 Operator: FRITSCHLAN\Benes

Material: wool ZrO2

P6, 15x20mm 60min

Description: water + 30s ultrasonic + 0,1% Na4P2O7

Calculation Automatische Modellerkennung

Mode Wet

Serial No. 22.2000.00/90771

Beam absorption 15,0 %

Pump 55,00 %

Meas. range 0,10 µm - 199,14 µm

Cellposition 2

Channels 102

Ultrasonic On

Scans 100

Prot.-No.: **090119**

5-99%

Obere Kornklasse [µm]	Q3(x) [%]
3,177	5,0
5,615	10,0
8,110	15,0
10,680	20,0
13,372	25,0
16,178	30,0
19,122	35,0
22,215	40,0
25,338	45,0
28,670	50,0
32,089	55,0
35,705	60,0
39,374	65,0
43,626	70,0
48,165	75,0
53,412	80,0
59,386	85,0
67,252	90,0
78,667	95,0
99,320	99,0

01-150µm

Obere Kornklasse [µm]	Q3(x) [%]
0,110	0,0
0,200	0,4
0,300	0,6
0,500	0,8
0,800	1,0
1,000	1,2
1,500	1,8
2,000	2,7
3,000	4,6
5,000	8,8
8,000	14,8
10,000	18,7
15,000	27,9
20,000	36,5
25,000	44,4
30,000	52,0
40,000	65,8
50,000	76,8
60,000	85,5
80,000	95,4
90,000	97,8
100,000	99,1
120,000	99,9
150,000	100,0

