

Meas. No. 9616 SOP 45 Date 12.05.2010 10:04:13 Operator: FRITSCHLAN\benes

Material: **Limestone**

**P5cl ZrO2 20mm 3min dry**

Description: 0,1% Na4P2O7 in water

Calculation Automatische Modellerkennung

Mode Wet

Serial No. 22.2000.00/90771

Beam absorption 13,6 %

Pump 60,00 %

Meas. range 0,10 µm - 331,90 µm

Cellposition 2

Channels 102

Ultrasonic On

Scans 100

Prot.-No.: **100135**

5-99%

Obere Kornklasse [µm]	Q3(x) [%]
0,644	5,0
0,910	10,0
1,188	15,0
1,504	20,0
1,885	25,0
2,361	30,0
2,992	35,0
3,799	40,0
4,822	45,0
6,050	50,0
7,478	55,0
9,111	60,0
11,040	65,0
13,435	70,0
16,561	75,0
20,984	80,0
27,208	85,0
35,724	90,0
47,735	95,0
67,354	99,0

01-100µm

Obere Kornklasse [µm]	Q3(x) [%]
0,110	0,0
0,200	0,1
0,300	0,4
0,500	2,6
0,800	7,9
1,000	11,7
1,500	19,9
2,000	26,4
3,000	35,1
5,000	45,8
8,000	56,7
10,000	62,5
15,000	72,7
20,000	79,0
25,000	83,5
30,000	86,8
40,000	92,0
50,000	95,7
60,000	98,0
80,000	99,8
90,000	99,9
100,000	100,0

