

Meas. No. 9624 SOP 45 Date 17.05.2010 10:44:07 Operator: FRITSCHLAN\benes

Material: **Limestone**

P5cl ZrO2 5mm 120min

Description: 0,1% Na4P2O7 in water

Calculation Automatische Modellerkennung

Mode Wet

Serial No. 22.2000.00/90771

Beam absorption 12,8 %

Pump 70,00 %

Meas. range 0,10 µm - 53,10 µm

Cellposition 1

Channels 51

Ultrasonic On

Scans 100

Prot.-No.: **100135**

5-99%

Obere Kornklasse [µm]	Q3(x) [%]
0,358	5,0
0,445	10,0
0,517	15,0
0,584	20,0
0,649	25,0
0,714	30,0
0,780	35,0
0,849	40,0
0,920	45,0
0,996	50,0
1,078	55,0
1,164	60,0
1,262	65,0
1,369	70,0
1,490	75,0
1,640	80,0
1,824	85,0
2,065	90,0
2,451	95,0
3,230	99,0

001-5 µm

Obere Kornklasse [µm]	Q3(x) [%]
0,010	0,3
0,050	0,3
0,100	0,3
0,200	0,6
0,300	2,5
0,400	7,2
0,500	13,7
0,600	21,2
0,700	28,9
0,800	36,5
0,900	43,6
1,000	50,2
1,500	75,4
2,000	88,8
2,500	95,4
3,000	98,3
3,500	99,5
4,000	99,9
4,500	100,0
5,000	100,0

