

Meas. No. 9639 SOP 45 Date 18.05.2010 09:23:04 Operator: FRITSCHLAN\benes

Material: **Limestone**

P5cl ZrO2 5mm 240min

Description: 0,1% Na4P2O7 in water

Calculation Automatische Modellerkennung

Mode Wet

Serial No. 22.2000.00/90771

Beam absorption 5,3 %

Pump 70,00 %

Meas. range 0,10 μm - 53,10 μm

Cellposition 1

Channels 51

Ultrasonic On

Scans 100

Prot.-No.: **100135**

5-99%

Obere Kornklasse [μm]	Q3(x) [%]
0,274	5,0
0,336	10,0
0,391	15,0
0,442	20,0
0,493	25,0
0,544	30,0
0,598	35,0
0,654	40,0
0,712	45,0
0,777	50,0
0,847	55,0
0,923	60,0
1,011	65,0
1,112	70,0
1,228	75,0
1,369	80,0
1,551	85,0
1,806	90,0
2,220	95,0
3,119	99,0

001-5 μm

Obere Kornklasse [μm]	Q3(x) [%]
0,010	0,1
0,050	0,1
0,100	0,1
0,200	1,1
0,300	7,0
0,400	15,9
0,500	25,7
0,600	35,2
0,700	44,0
0,800	51,8
0,900	58,5
1,000	64,4
1,500	83,8
2,000	92,7
2,500	96,9
3,000	98,7
3,500	99,5
4,000	99,8
4,500	99,9
5,000	100,0

