

Meas. No. 9622 SOP 45 Date 12.05.2010 15:38:28 Operator: FRITSCHLAN\benes

Material: **Limestone**

**P5cl ZrO2 5mm 60min**

Description: 0,1% Na4P2O7 in water

Calculation Automatische Modellerkennung

Mode Wet

Serial No. 22.2000.00/90771

Beam absorption 13,2 %

Pump 70,00 %

Meas. range 0,10 µm - 53,10 µm

Cellposition 1

Channels 51

Ultrasonic On

Scans 100

Prot.-No.: **100135**

5-99%

Obere Kornklasse [µm]	Q3(x) [%]
0,443	5,0
0,558	10,0
0,658	15,0
0,753	20,0
0,849	25,0
0,946	30,0
1,047	35,0
1,153	40,0
1,266	45,0
1,388	50,0
1,519	55,0
1,665	60,0
1,831	65,0
2,015	70,0
2,232	75,0
2,507	80,0
2,840	85,0
3,300	90,0
4,026	95,0
5,627	99,0

001-10µm

Obere Kornklasse [µm]	Q3(x) [%]
0,015	0,0
0,040	0,0
0,050	0,0
0,100	0,0
0,200	0,0
0,300	0,9
0,500	7,3
0,800	22,4
1,000	32,7
1,500	54,4
2,000	69,6
3,000	86,9
4,000	94,9
5,000	98,1
6,000	99,3
7,000	99,8
8,000	100,0
9,000	100,0
10,000	100,0

