

**Meas.Nr.** 71078      **Date** Donnerstag, 05 Dezember 2013 08:27      **Personal Report No.:**  
**Material** Asbestos      60 min P7pl 3 mm ZrO2 in water      **M130287**  
**Description** in water + 0,1 % Na4P2O7 - 1 min Ultrasonic LABORETTE 17

**Calculation** Fraunhofer      broad      **TradeOff** 1.000000e+002  
**Scans Fine** 100      **Scans Coarse** 0      **Channels** 102      **Beam Obscuration** 10.00 %  
**Meas. Range** 0.01 [µm] - 42.30 [µm]      **Pump** 3      **Ultrasonics** 10

%	< µm
5	0.25
10	0.31
15	0.38
20	0.44
25	0.51
30	0.59
40	0.77
50	1.00
60	1.29
70	1.67
80	2.20
90	3.13
95	4.17
98	6.92
99	10.17

%	< µm
0.1	0.10
24.0	0.50
50.0	1.00
65.9	1.50
76.6	2.00
89.0	3.00
94.5	4.00
96.7	5.00
97.6	6.00
98.0	7.00
98.4	8.00
99.0	10.00
99.4	12.00
99.9	15.00
100.0	20.00

