

Meas. No. 7817 SOP 36 Date 08.10.2009 09:09:46 Operator FRITSCHLAN\Benes

Material: BT

P7pl 80ml ZrO2 5mm 20min

Description: dispersed in isopropyl alcohol (2min ultrasonic)

Calculation Automatische Modellerkennung

Mode SVDU

Serial No. 22.2000.00/90771

Beam absorption 13,0 %

Pump ca. 60 %

Meas. range 0,10 µm - 53,10 µm

Cellposition 1

Channels 51

Ultrasonic Off

Scans 100

Prot.-No.: 090243

5-99%

Obere Kornklasse [µm]	Q3(x) [%]
0,567	5,0
0,697	10,0
0,805	15,0
0,905	20,0
1,002	25,0
1,101	30,0
1,202	35,0
1,307	40,0
1,421	45,0
1,543	50,0
1,676	55,0
1,831	60,0
2,003	65,0
2,206	70,0
2,474	75,0
2,819	80,0
3,347	85,0
4,389	90,0
7,105	95,0
11,350	99,0

01-20µm

Obere Kornklasse [µm]	Q3(x) [%]
0,100	0,1
0,150	0,4
0,200	0,6
0,300	1,0
0,400	1,5
0,500	3,1
0,600	6,1
0,700	10,1
0,800	14,8
1,000	24,9
1,200	34,9
1,500	48,4
2,000	64,9
2,500	75,4
3,000	81,9
4,000	88,7
5,000	91,6
6,000	93,4
8,000	96,2
10,000	98,2
15,000	99,9
20,000	100,0

