

Meas. No. 7815 SOP 36 Date 08.10.2009 08:52:11 Operator FRITSCHLAN\Benes

Material: LZZ

P7pl 80ml ZrO2 5mm 20min

Description: dispersed in isopropyl alcohol (2min ultrasonic)

Calculation Automatische Modellerkennung

Mode SVDU

Serial No. 22.2000.00/90771

Beam absorption 11,7 %

Pump ca. 60 %

Meas. range 0,10 µm - 185,86 µm

Cellposition 2

Channels 102

Ultrasonic Off

Scans 100

Prot.-No.: 090243

5-99%

Obere Kornklasse [µm]	Q3(x) [%]
0,124	5,0
0,154	10,0
0,187	15,0
0,226	20,0
0,270	25,0
0,319	30,0
0,379	35,0
0,452	40,0
0,541	45,0
0,651	50,0
0,787	55,0
0,952	60,0
1,145	65,0
1,373	70,0
1,641	75,0
1,971	80,0
2,389	85,0
3,024	90,0
4,250	95,0
8,413	99,0

01-20µm

Obere Kornklasse [µm]	Q3(x) [%]
0,100	1,7
0,150	9,4
0,200	16,7
0,300	28,2
0,400	36,6
0,500	42,8
0,600	47,8
0,700	51,9
0,800	55,4
1,000	61,3
1,200	66,3
1,500	72,5
2,000	80,4
2,500	86,0
3,000	89,9
4,000	94,4
5,000	96,5
6,000	97,7
8,000	98,8
10,000	99,5
15,000	100,0
20,000	100,0

