

Meas. No. 1975    SOP 55    Date 06.11.2007 15:59:33    Operator: FRITSCHLAN\benes

Material: **Mineral**

**60min premium line**

Description: 60min premium line 80ml ZrO2 with 100g 1-3mm balls

Calculation Automatische Modellerkennung

Mode Wet

Serial No. 22.2000.00/90771

Beam absorption 10.0 %

Pump 40.00 %

Meas. range 0.01  $\mu\text{m}$  - 53.08  $\mu\text{m}$

Cellposition 1

Channels 108

Ultrasonic An

Scans 100

Prot.-No.:

5-99%

Obere Kornklasse [ $\mu\text{m}$ ]	Q3(x) [%]
0.169	5.0
0.188	10.0
0.203	15.0
0.216	20.0
0.229	25.0
0.240	30.0
0.252	35.0
0.264	40.0
0.276	45.0
0.289	50.0
0.303	55.0
0.316	60.0
0.333	65.0
0.350	70.0
0.370	75.0
0.391	80.0
0.424	85.0
0.468	90.0
0.549	95.0
1.547	99.0

001-5  $\mu\text{m}$

Obere Kornklasse [ $\mu\text{m}$ ]	Q3(x) [%]
0.010	0.0
0.050	0.0
0.100	0.0
0.200	13.8
0.300	54.1
0.400	81.5
0.500	92.7
0.600	96.2
0.700	97.0
0.800	97.1
0.900	97.2
1.000	97.3
1.500	98.9
2.000	99.8
2.500	100.0
3.000	100.0
3.500	100.0
4.000	100.0
4.500	100.0
5.000	100.0

